



## **Cirrus HD-OCT**

Certainty in seconds. Certainty for years.



We make it visible.

// CERTAINTY  
MADE BY CARL ZEISS

## We know you'll love a Cirrus.

Keeping both your patients and your practice in mind, Carl Zeiss Meditec, the global leader in OCT, developed Cirrus™ HD-OCT. Not only does it supply you with bar-setting imagery, it delivers detailed diagnostic and change analyses you can rely on time and again. Along with its small footprint and fast capture speeds, Cirrus is designed to improve workflow efficiency while helping you deliver better care to your patients.

*It's time to see what you've been missing. It's time for Cirrus.*

### Superior analysis

With high-density cube data and proven segmentation, Cirrus delivers a diagnostic analysis you can trust.

#### Central Subfield Macular Thickness Repeatability Standard Deviation

No Disease	2.5 $\mu$ m
AMD	8.7 $\mu$ m
Macular Edema	7.0 $\mu$ m
Diabetic Retinopathy	8.1 $\mu$ m
VRI Disorder	4.3 $\mu$ m

Source: 510(k) Summary: Cirrus HD-OCT with Retinal Nerve Fiber Layer and Macular Normative Databases, [www.accessdata.fda.gov/cdrh\\_docs/pdf8/K083291.pdf](http://www.accessdata.fda.gov/cdrh_docs/pdf8/K083291.pdf).

#### Receiver Operating Characteristic Curves, Normal vs. Glaucomatous Eyes


Parameter	Mild	Moderate to Severe
Average RNFL Thickness	0.893	0.993
Rim Area	0.912	0.999
Vertical Cup-to-Disc Ratio	0.890	0.995

Source: Mwanza et al. Ability of Cirrus HD-OCT optic nerve head parameters to discriminate normal from glaucomatous eyes. *Ophthalmology*. 2011;118(2):241-248.

### Spectacular imagery

With legendary ZEISS optics and Cirrus, you'll experience brilliant, detail-rich visuals to help you diagnose and care for your patient.


Cirrus uses Selective Pixel Profiling™ to optimize each pixel in its HD Raster Scans. It produces imagery that goes beyond mere image-averaging. It's a difference you need to see to believe.




## One pathology. Multiple views.

Cirrus gives you the ability to view pathologies from multiple vantage points—and with a range of at-a-glance visualization formats, you'll be able to better assess your patient's condition and determine the appropriate course of action.


## Postsurgical Pseudophakic Cystoid Macular Edema



HD 5 Line Raster




LSO Fundus Image with Raster Line



Macular Thickness Map


## Age-Related Macular Degeneration



HD 5 Line Raster




RPE Segmentation Map




Macular Thickness Map


## Epiretinal Membrane




HD 5 Line Raster




ILM Segmentation Map




Macular Thickness Map




LSO Fundus Image with Macular Thickness Map



Advanced Visualization™ with RPE Fit Slab




3D Visualization



LSO Fundus Image with ILM Slab




Advanced Visualization™ with ILM Slab



3D Visualization



## Discover the power of the Cirrus Cube.



Cirrus offers unsurpassed OCT technology. Capturing a tightly packed, detail-rich cube of data in just seconds, it allows you to both visualize and analyze your patient's condition. Because the cube is populated with such high-density data, you can explore pathologies without requiring additional scan patterns.

Scan Pattern	Data Points Per A-Scan	Total Data Points	Spacing Between Lines	Capture Time
512 x 128	1024	> 67 million	47 $\mu\text{m}$	2.4 s
200 x 200	1024	> 40 million	30 $\mu\text{m}$	1.5 s

### Scan with greater granularity

Closely spaced B-scans within the cube ensure that even small areas of pathology are captured and easily viewable, unlike scans that are spaced further apart, which may miss the central fovea or nearby subtle defects.

### Enhance your analysis

Millions of data points from the cube are fed into ZEISS proprietary algorithms for accurate segmentation, reproducible measurements and registration for change analysis.

## Analysis you can trust.

Generating a comprehensive cube of data is only the beginning. Cirrus gives you the ability to see beyond the scan and transform information into insight, becoming an indispensable part of your day-to-day clinical decision-making process.

### Algorithm excellence

Carl Zeiss Meditec and its research collaborators have developed advanced algorithms to measure and display layers.

### Cube registration to track change

Cirrus data cubes are automatically registered with data from prior visits, allowing for point-to-point comparisons.

### Automatic, accurate centering of the measurement

FoveaFinder™ and AutoCenter™ technologies ensure that measurements are made in the correct locations, taking the pressure off the operator to center the scans perfectly.



### Normative data

Diversified normative databases for ONH, RNFL and macular thickness facilitate at-a-glance assessments.



### Increase your diagnostic certainty

Cirrus enables rapid, careful assessment of the retina. By utilizing precise macular thickness analyses, providing detailed ILM and RPE layer maps and putting more than 100 B-scans at your disposal, Cirrus provides the framework to assess your patient's retinal condition.




#### FoveaFinder™

With FoveaFinder™, Cirrus automatically and accurately locates the fovea, centers the ETDRS measurement grid and presents you with the B-scan through the fovea by default.

### Track subtle macular change

Cirrus data cubes are automatically registered with data from prior visits after the scan is acquired. This enables side-by-side visualization of the same location on the retina for each visit. Cirrus compares measurements from the current and prior visits to provide a thickness change map and helps you determine next steps for your patient.






# // GLAUCOMA MANAGEMENT

MADE BY CARL ZEISS

## Identify and track RNFL and ONH for glaucoma management

With Cirrus, all traditional RNFL measurements based on the 3.4 mm circle are present; however, Cirrus enables you to see past the circle-based assessments. Spotting wedge defects and other patterns of loss is simplified with Deviation Maps, which show comparisons to normative data for each superpixel in the 6 x 6 mm area.

Unique Cirrus Optic Nerve Head analysis provides automated identification of the optic disc and cup boundaries. The analysis is generated using the dense data in the Optic Disc 200 x 200 data cube in tandem with a proprietary ZEISS algorithm. This algorithm precisely measures the neuroretinal rim while accounting for tilted discs, disruptions to the RPE and other challenging pathologies.




### AutoCenter™


After the scan is acquired, Cirrus automatically centers the measurement circle around the disc. The placement is not operator-dependent.



**Guided Progression Analysis™ (GPA™)** compares RNFL thickness measurements from data cubes obtained during different visits and allows you to determine if statistically significant change has occurred over time.



**The HFA-Cirrus Combined Report**, available exclusively with ZEISS FORUM®, summarizes patient structure and function information in a single display.




# // ANTERIOR SEGMENT IMAGING

MADE BY CARL ZEISS




## Expand your diagnostic insight


Cirrus offers anterior segment imaging of the angle and cornea and the ability to measure central cornea thickness.




*Narrow Angle Visualization*



*Central Corneal Thickness Measurement*




*Keratic Precipitates*



*Microstriae*




*PRK Scar*




*Descemet Detachment*



*IEK Zigzag*



*Bullous Keratopathy with BCL*



*Cataract Incision*

**From the industry leader in OCT, Cirrus is the best-selling spectral domain OCT in the world.**

Cirrus represents the culmination of decades of patents, prototypes and progress. ZEISS is committed to delivering the excellence in installation, training and ongoing support you expect from the market leader.

As new diagnostic needs emerge and new therapies are developed, innovation continues with Cirrus.

In addition, recognizing the modern electronic workplace, Cirrus integrates seamlessly into EMRs and with FORUM<sup>®</sup>, our advanced data management solution for simplifying, centralizing and viewing the vast amounts of clinical data generated by ophthalmic instruments.

There's only one OCT that promises you Certainty in Seconds, Certainty for Years.<sup>™</sup>

There's only one Cirrus.



For videos, presentations, recent clinical literature and updated product information, visit:

[www.meditec.zeiss.com/cirrus](http://www.meditec.zeiss.com/cirrus)

**Your local contact:****Argentina**

Carl Zeiss Argentina S.A.  
Calle Nahuel Huapi 4015 / 25  
C1430 BCO Buenos Aires  
Argentina  
Phone: +54 11 45 45 66 61  
bruzzi@zeiss.com.ar

**Australia**

Carl Zeiss Pty. Ltd.  
Unit 13, 2 Eden Park Drive  
North Ryde, New South Wales  
2113  
Australia  
Phone: +61 2 9020 1333  
med@zeiss.com.au

**Austria**

Carl Zeiss GmbH  
Laxenburger Str. 2  
1100 Vienna  
Austria  
Phone: +43 1 79 51 80  
austria@zeiss.org

**Belgium**

Carl Zeiss NV-SA  
Ikaroslaan 49  
1930 Zaventem  
Belgium  
Phone: +32 2 719 39 11  
info@zeiss.be

**Brazil**

Carl Zeiss do Brasil Ltda.  
Av. Nações Unidas, 21711  
CEP04795-100 São Paulo  
Brazil  
Phone: +55 11 5693 5521  
medbrasil@zeiss.org

**Canada**

Carl Zeiss Canada Ltd.  
45 Valleybrook Drive  
Toronto, ON M3B 2S6  
Canada  
Phone: +1 800 387 8037  
micro@zeiss.com

**China**

Carl Zeiss Shanghai Co. Ltd.  
1/f., Ke Yuan Building  
11 Ri Yin Nan Road  
Waigaoqiao Free Trade Zone  
2005 Yang Gao Bei Road  
Shanghai 200131  
China  
Phone: +86 21 5048 17 17  
sro@zeiss.com.cn

**Czech Republic**

Carl Zeiss spol. s.r.o.  
Radlická 14/3201  
150 00 Prague 5  
Czech Republic  
Phone: +420 233 101 221  
zeiss@zeiss.cz

**France**

Carl Zeiss Meditec France SAS  
60, route de Sartrouville  
78230 Le Pecq  
France  
Phone: +33 1 34 80 21 00  
med@zeiss.fr

**Germany**

Carl Zeiss Meditec VG mbH  
Carl-Zeiss-Strasse 22  
73446 Oberkochen  
Germany  
Phone: +49 7364 20 6000  
vertrieb@meditec.zeiss.com  
Surgical Ophthalmology:  
Phone: +49 800 470 50 30  
iol.order@meditec.zeiss.com

**Hong Kong**

Carl Zeiss Far East Co. Ltd.  
Units 11-12, 25/F  
Tower 2, Ever Gain Plaza  
No. 88 Container Port Road  
Kwai Chung  
Hong Kong  
Phone: +852 2332 0402  
czfe@zeiss.com.hk

**India**

Carl Zeiss India Pvt. Ltd.  
22, Kensington Road  
Ulsoor  
Bangalore 560 008  
India  
Phone: +91 80 2557 88 88  
info@zeiss.co.in

**Italy**

Carl Zeiss S.p.A.  
Viale delle Industrie 20  
20020 Arese (Milan)  
Italy  
Phone: +39 02 93773 1  
post@zeiss.it

**Japan**

Carl Zeiss Meditec Japan Co. Ltd.  
Shinjuku Ku  
Tokyo 160-0003  
22 Honchio-Cho  
Japan  
Ophthalmic instruments:  
Phone: +81 3 33 55 0331  
medsales@zeiss.co.jp  
Surgical instruments:  
Phone: +81 3 33 55 0341  
cmkoho@zeiss.co.jp

**Malaysia**

Carl Zeiss Sdn Bhd.  
Lot 2, Jalan 243/51 A  
46100 Petaling Jaya  
Selangor Darul Ehsan  
Malaysia  
Phone: +60 3 7877 50 58  
malaysia@zeiss.com.sg

**Mexico**

Carl Zeiss de México S.A. de C.V.  
Avenida Miguel Angel de Quevedo 496  
04010 Mexico City  
Mexico  
Phone: +52 55 59 99 0200  
cz-mexico@zeiss.org

**Netherlands**

Carl Zeiss B.V.  
Trapezium 300  
Postbus 310  
3364 DL Sliedrecht  
Netherlands  
Phone: +31 184 43 34 00  
info@zeiss.nl

**New Zealand**

Carl Zeiss (N.Z.) Ltd.  
15B Paramount Drive  
P.O. Box 121 - 1001  
Henderson, Auckland 0650  
New Zealand  
Phone: +64 9 838 5626  
med@zeiss.com.a

**Poland**

Carl Zeiss sp. z o.o.  
ul. Lopuszanska 32  
02-220 Warsaw  
Poland  
Phone: +48 22 858 2343  
medycyna@zeiss.pl

**Singapore**

Carl Zeiss Ptd. Ltd.  
50 Kaki Bukit Place  
Singapore 415926  
Singapore  
Phone: +65 6741 9600  
info@zeiss.com.sg

**South Africa**

Carl Zeiss (Pty.) Ltd.  
363 Oak Avenue  
Ferndale  
Randburg 2194  
South Africa  
Phone: +27 11 886 9510  
info@zeiss.co.za

**South Korea**

Carl Zeiss Co. Ltd.  
Seoul 121-828  
Mapo-gu  
141-1, Sangsu-dong  
2F, BR Elitel Bldg.  
South Korea  
Phone: +82 2 3140 2600  
korea@zeiss.co.kr

**Spain**

Carl Zeiss Meditec Iberia S.A.  
Ronda de Poniente, 15  
Tres Cantos  
28760 Madrid  
Spain  
Phone: +34 91 203 37 00  
info@zeiss.es

**Sweden**

Carl Zeiss AB  
Tegeluddsvaegen 76  
10254 Stockholm  
Sweden  
Phone: +46 84 59 25 00  
info@zeiss.se

**Switzerland**

Carl Zeiss AG  
Feldbachstrasse 81  
8714 Feldbach  
Switzerland  
Phone: +41 55 254 7534  
med@zeiss.ch

**Thailand**

Carl Zeiss Thailand  
Floor 8, Thosapol Land Building 2  
230 Ratchadapisek Road  
Huaykwang, Bangkok 10310  
Thailand  
Phone: +66 2 2 74 06 43  
thailand@zeiss.com.sg

**United Kingdom**

Carl Zeiss Ltd.  
15-20 Woodfield Road  
Welwyn Garden City  
Hertfordshire, AL7 1JQ  
United Kingdom  
Phone: +44 1707 871200  
info@zeiss.co.uk

**United States of America**

Carl Zeiss Meditec, Inc.  
5160 Hacienda Drive  
Dublin, CA 94568  
USA  
Phone: +1 925 557 4100  
info@meditec.zeiss.com



**Carl Zeiss Meditec, Inc.**

5160 Hacienda Drive  
Dublin, CA 94568  
USA  
www.meditec.zeiss.com/cirrus



**Carl Zeiss Meditec AG**

Goeschwitzer Str. 51-52  
07745 Jena  
Germany  
www.meditec.zeiss.com/cirrus